




Cardiology Specialists of Acadiana

315 Rue Louis XIV, Lafayette, LA 70508

P 337-269-9777 • F 337-269-0244


Patient Instructions: Holter & Cardiac Monitor

24 Hour



5 Lead Electrode Placement

14 Day



24 Hour Patient Instructions:

- If a cable comes off, snap it back on.
- Do not remove stickers.
- Do not tamper with the electrodes, wires, or recorder.
- Do not use an Electric Blanket.
- DO NOT GET MONITOR WET.
- Make sure you hit the event button if you have symptoms.
- Do not remove recorder from its case.

14 Day Patient instructions:

- You may remove monitor to take bath then immediately put back on.
- Do not remove stickers unless loose.
- Do not pull on wires to remove leads. Grab by the end
- If sticker falls off replace with a new one (provided)
- Follow diagram below to snap leads back on
- DO NOT GET MONITOR WET
- Make sure you hit the event button if you have symptoms
- When reapplying electrodes try to place them on the collar bone and on a rib

